


Laryngeal Mask Airway

Latex
Free



Reusable Silicone LMA



Disposable Silicone LMA



Disposable PVC LMA

- Used for **Inhalational anaesthesia**
- **Emergency management** of airway in failed intubation
- Perfect for **supporting** and **maintaining** the airway without tracheal intubation
- Allows **single-handed** ventilation
- Allows **rapid & blind** insertion (no laryngoscopy required)
- **No** muscle relaxants required
- A **secure, hands free** airway
- **Smooth emergence** from anaesthesia
- **Minimal** soft tissue and dental **trauma**
- **Latex-free**
- Available in **Silicone** or **PVC**

Size	Product Code		
	Reusable Silicone	Disposable Silicone	Disposable PVC
1.0	509-008-70	509-008-10	509-009-10
1.5	509-008-75	509-008-15	509-009-15
2.0	509-008-80	509-008-20	509-009-20
2.5	509-008-85	509-008-25	509-009-25
3.0	509-008-90	509-008-30	509-009-30