

Ulma INDICATOR

Steam Autoclave Tape

Features:

- It provides an economical method of determining sterility
- The versatile tape adheres to all wraps, packs and sets. It can be used with most type of materials – paper, plastic, woven, non-woven
- It doubles up as a sealing tape for wraps and packs being sterilised. The strong tape holds firmly to packs before, during and after sterilisation and can be removed easily without leaving any residue
- The paper surface of tape allows writing on it for good record keeping or making notes

Applications:

- Monitor sterilisation exposure in steam sterilisation process
- The chemical process indicator lines on the tape change from white to black on exposure to steam to provide visual evidence of sterilisation
- Secure wraps, packs or sets containing items which need to be sterilised
- Tape can be applied to various locations of the packs to ensure even and complete sterilisation



white lines change to black

Size	Packing (rolls/box)	Product Number
12mm (1/2") x 55m/roll	12	516-001-10
19mm (3/4") x 55m/roll	8	516-001-20
25mm (1") x 55m/roll	6	516-001-30

* Custom widths and lengths are available upon request