

Accoson Greenlight 300

accuracy confidence

self-calibration without mercury

Self-calibration to zero each time the device is switched on ensures reliable accuracy

The strength of the greenlight 300 is that it is suitable for any patient group (including the pregnant and the elderly) where mercury was used in the past

This is because:

- It is a manual device
- It displays every 2mmHg
- It is easier to read than mercury
- It is calibrated to 0.4mmHg
- It carries out an auto calibration to zero each time it is switched on
- The cuff deflation rate is displayed
- Any cuff size can be used
- Frequent calibration checks are not required. There is none of the fragility of an aneroid gauge

Wall model The greenlight 300 converts to a wall mounted instrument when fitted with the Wall Conversion Kit Part Number 505-018-06 consisting of:



- A swivel mounting bracket which allows the instrument to be turned in two directions. This gives easy viewing of the display particularly when the instrument is placed between two beds
- A separate cuff container to accommodate the cuff and tubing

Stand Model the greenlight 300 converts to a mobile stand model when fitted with the Stand Conversion Kit 505-018-07 consisting of:

- A large container for the cuff and tubing on which the greenlight 300 is fitted. The instrument is set at an angle for easy viewing
- A mounting stem and a heavy cast base with five castors to give mobility and stability



Greenlight 300 with Adult Velcro cuff	505-018-01
Greenlight 300 with Wide Range Adult Velcro cuff	505-018-02
Greenlight 300 Gauge only	505-018-03
Greenlight 300 with Stethoscope combined in cuff	505-018-04
Greenlight 300 Calibration Test Kit	505-018-05
Wall conversion kit	505-018-06
Stand conversion kit	505-018-07
Carrying case with zipper	505-018-08
Set of 4 x AA Batteries	505-018-09