

Hydrocolloid Dressing



Description:

- Consist of non adhesive Hydrocolloid padding with transparent film backing for monitoring

Features:

- Easy to apply, does not cause trauma on removal
- Hypoallergic
- Transparent film backing, enables monitoring of fluid absorption for timely dressing changes
- Dressing ensures an optimal moist environment, comfortable for extended wear time.
- Available in Regular and Thin dressings

Application:

- The Regular Hydrocolloid Dressing is effective for management of light to moderately exuding wounds as well as first and second degree burns
- The Thin Hydrocolloid Dressing is suitable for management of light exuding wounds

Specification:

- Sterile
- With natural latex
- Available in oval and square shapes

Packing:

Size	Pad Size	Shape	Pcs/Box	Product Code
10cm x 12cm	7cmx 9cm	Oval	10pcs/box	500-020-01
13cm x 15cm	10cmx12cm	Oval	10pcs/box	500-020-02
17cm x 20cm	14cmx17cm	Oval	6pcs/box	500-020-03
10cm x 10cm	10cmx10cm	Square	5pcs/box	500-020-04
15cm x 15cm	15cmx15cm	Square	5pcs/box	500-020-05
20cm x 20cm	20cmx20cm	Square	5pcs/box	500-020-06