


Mesh Dressing with Silver



Description:

- For fast acting, long lasting, barrier control.
- Provides an effective antimicrobial barrier in silver against a wide range of microbes including antibiotic-resistant bacteria such as MRSA and VRE

Features:

- Soft, comfortable, ready to use
- No minimal staining and non irritating
- Contains silver sulfate which is released as silver ions in the dressing and creates an effective antimicrobial barrier for up to seven days
- The porous, non-occlusive dressing conforms to the wound base and wicks drainage into the dressing where the silver ions are available to reduce the number of bacteria and yeast

Application:

- For use as a wound dressing over abrasions, ulcers, trauma wounds, surgical wounds, first and second degree burns and donor sites.
- Ideal for use in wound clinics, hospitals, home care and long term care

Specification:

- Non Woven, porous cotton
- Contains silver sulfate (less than 1%)
- Soft and highly absorbent.
- Sterile

Packing:

Size	Pcs/Box	Box/Carton	Product Code
5cm x 5cm	5 pcs/box	4 boxes/carton	500-022-01
10cm x 12cm	5 pcs/box	4 boxes/carton	500-022-02
10cm x 20cm	5 pcs/box	4 boxes/carton	500-022-03
20cm x 20cm	3 pcs/box	4 boxes/carton	500-022-04
40cm x 40cm	3 pcs/box	4 boxes/carton	500-022-05